

● MAGNETIC STIRRER

Design to mix watery liquid, provided with soft Start facility, electronic speed regulator & fitted with S.S. Top and supplied Teflon Magnetic Bit.



Types of MAGNETIC STIRRERS:

- » Magnetic Stirrer
- » Multi Magnetic Stirrers
(Available with 2 in 1, 3 in 1, 6 in 1)
- » Magnetic Stirrer with Hot Plate
- » Multi Magnetic Stirrers with Hot Plates
(Available with 2 in 1, 3 in 1, 6 in 1)

Salient Features

- » Heavy duty PMDC Motor & Magnet.
- » Heavy duty solid state electronic step less Speed Regulator.
- » Max. Temperature up to 200°C for Magnetic Stirrer with HOTPLATE

Optional Accessories

- » Digital Rpm Indicator / Controller
- » Digital Temperature Indicator cum Controller for Hot Plate Model
- » Digital Switch off Timer

Liquid stirring Capacity

0.5 to 1 Liter

2 Liter


5 Liter



with DIGITAL RPM INDICATOR

DIGITAL RPM INDICATOR & TEMPERATURE CONTROL

● MAGNETIC STIRRER WITH HOT PLATE



MANUAL CONTROL



DIGITAL TEMPERATURE CONTROLLER